

Manufacturer	Medical Device	Model
Axia	Patient Monitor	V1050T, V1200T, VT1500A, i120
Asepct Medical Systems (Covidien)	Brain Monitoring System	A-2000 BIS, A-2000 XP BIS, VISTA BIS
BioLight	Patient Monitor	Q2, Q3, Q4, Q5, Q6, Q7
Blink	Neuromuscular Monitor	TwitchView
Criticare Systems	Patient Monitor	8100H nCompass, 8500H Poet IQ, 8500Q Poet IQ2 506 eQuality Series, 505 VitalCare Series, nGenunity Series
CAS Medical	Patient Monitor	740 ,740 Select, 750 FORE-Sight, FORE-Sight Elite
Drager Medial	Patient Monitor	Gamma, Gamma XL, Gamma X XL Delta, Delta XL, Kappa, Kappa XLT Vista, Vista XL, Vista 120, Vista 120S, Infinity ACS
	Anesthesia Machine	Apollo, Babylog (8000, 8000 SC, 8000 Plus, VN500) Carina, Cato, Cicero B, Cicero C, Cicero EM Evita, Evita 2, Evita 4, Evita Infinity V500, Evita V300, Evita XL Fabius (GS, OS, MRI, Plus, Trio), Julian, Pallas, Perseus A500 Primus, PhysioFlex, SA2, Savina, Sulla 808 V, Titus, Trajan Vamos, Zeus, Atlan (A300, A350)
DRE	Patient Monitor	Waveline Pro (V1), Waveline Pro (V2) Waveline EZ, Waveline EZ MAX
Edan	Patient Monitor	Elite V5, V6, V8 iM50, iM60, iM70, iM80, M50, M80, x8, x10, x12
Edward Lifesciences	Hemodynamic Monitor	EV1000, Vigilance, Vigilance II, Vigileo, HemoSphere
GE Healthcare (Datex Ohmeda)	Patient Monitor	B30 ¹ , B40 ¹ , B105 ¹ , B125 ¹ , B105M, B125M Cardiacap II, Cardiacap/5 Carescape B450, B650, B850, V100 Dash 2000, 2500, 3000, 4000, 5000 Anesthesia Monitor AS/3, S/5 Critical Care Monitor CS/3, S/5 FM Monitor S/5, FM Light Monitor S/5 Solar 7000, 8000, 8000i, 8000M, 9500 Procure400
	Anesthesia Machine	Aestiva (3000, 7100, 7900, Aestiva/5, 7800, MRI) Non-Aestiva 7900, Aespire (7900, View) Aisys, Aisys CS2, Avance, Avance CS2 Carestation (620, 650, 750, 850) Centiva /5, ADU-98, Engstrom Carestation. iVent 201 Modulus CD

¹ Requires the GE Multi I/O adapter, please refer to the supplemental note

Manufacturer	Medical Device	Model
IDMed	Neuromuscular Monitor	ToFscan
Invivo	MRI Monitor	Expression
Iradimed	MRI Monitor	3880 MRI Monitor
Larsen & Toubro	Patient Monitor	Star 55
Masimo	Patient Monitor	Radical, Root
Maquet/Getinge	Anesthesia Machine	Flow-i
Medred	MRI Monitor	Veris
Medtronic	Patient Monitor	Invos, Capnostream 20p
	Brain Monitoring System	BIS Advance
Mindray	Patient Monitor	Expert, MDPro Guardian Plus Passport 2 ² , 8, 12, 17m, XG, V Accutorr Plus, 3, 7, V Spectrum, Spectrum OR Trio (only for SN 0988-00-0600-4xxx) Beneview T1, T5, T8, N22, N17, N15, N12, N1 VS8, VS9, V12, V21, DPM6, DPM7, cPM8, cPM12 ePM12M, ePM12, ePM10, ePM10A, ePM10M, ePM15M
	Anesthesia Machine	A3, A4, A5, A7, A8, A9, AS3000
Mortara	Patient Monitor	SurveyorS12, SurveyorS19
Nihon Kohden	Patient Monitor	BSM-2301, 2304, 2351, 2353 BSM-3104 (with pre-installed QI-373P in the monitor) BSM-5105, 5106, 5135, 5136 BSM-6301, 6311, 6312, 6313, 6331, 6332, 6333, 6501, 6701 ³ BSM-9172 TR-6000 Series, TR-9000 Series, Life Scope G5, G5 Max, G7, G9
	Patient Monitor	K, Omni II
Omni	Patient Monitor	K, Omni II
Paragon	Anesthesia Machine	Platinum SC430
Perimed	Patient Monitor	P5000
Philips Healthcare	Patient Monitor	CMS, MP2, MP5 Intellivue ⁴ (MP20-MP90, MX450-MX850), MMS X2 SureSigns VS1, VS2+, VS3, VS4, VSi, VM1, MV4, MV6, VM8 Expression IP5, MR400, MR Patient Care Portal 5000
	Patient Monitor	TCM TOSCA
Radio Meter	Patient Monitor	TCM TOSCA
Skarray	Patient Monitor	Star65, Star90
SpaceLabs Healthcare	Patient Monitor	Qube, Xprezzon Ultraview 1500, 1600, 1700 Ultraview SL 2200, 2400, 2600, 2700, 2800, 2900, 3800, 3900
	Anesthesia Machine	Blease Sirius

² Requires an open SP1 or SP2 port, please refer to the supplemental note

³ Requires an QI-61P board installed , please refer to the supplemental note


⁴ Intellivue monitors must have an available LAN, MIB or RS232 port

Supplemental Note

GE B20, B40

Requires the GE Multi I/O adapter. GE has discontinued this adapter.


B40 and B20 Multi IO Adapter P/N: 2053489-051



GE B105, B125

Requires the GE Multi I/O adapter. Please purchase it from a GE.



B125 and B105 Multi IO Adapter P/N: 2098800-018



Mindray Passport 2

Please make sure there is a port call SP1 or SP2 and it is not used by a gas module.

RD1 is the same 9-pin male. But it does not do data export.



Nihon Kohden 67xx

Please make sure there is a DB9 port (9 pins, 2 rows). If the port is missing you will need to purchase a QI-671P board from Nihon Kohden

